

- EN** Alternating Therapy kit mattress and pump
F Kit matelas anti-escarres et compresseur à cycle alterné
ES Kit colchón antiescaras y compresor de aire alternante
P Kit colchão anti-escaras e compressor alterno

píuma
up2

Interchangeable cells mattress and pump with pressure regulation

Cod. **LAD463** *Kit interchangeable cells mattress and pump with pressure regulation*

Cod. **LAR053** *Only mattress in Nylon + PVC*

Cod. **LAR243** *Only Piuma Up compressor with regulation*

FUNCTIONS



ALTERNATING



Suitable for the prevention and treatment of **II STAGE** decubitus sores. Alternating cycle 12 min. Pump with pressure regulation according to the patient's weight. Cells fixed with buttons, easy to change in case of damage, and washable. PU or stretch cover **OPTIONAL**.

Technical features

SPECIFICATION	
Mattress	
Material	Nylon®/PVC
Material thickness	0,3 mm
Cell height	About 10,2 cm / 4"
Number of cells	20
Flaps	50 cm
Compressor	
Power	220V / 50-60 Hz
Pressure	40-100 mmHg
Power cable	2 m
Air release	6-7 l/mn
Cycle time	12 minutes
Noise	32 - 34 dB

WEIGHT AND DIMENSION					
	W:	D:	H:	Kg	Icon
Compressor	27 cm	13 cm	10,5 cm	1 kg	-
Mattress without flaps					
While functioning	200 cm	86 cm	10,2 cm	3,8 kg	145 kg
Spare parts			Optional		
 Cod. LAR204 Spare cell LAR053			 Cod. LAR304 PU cover for LAR053		
			 Cod. LAR305 Stretch cover for LAR053		